




















Disposable Air/Water, Suction, Biopsy Valve & Auxiliary Water Connector Kit.

Disposable Air/Water Channel Cleaning Adapter and Air/Water, Suction, Biopsy Valve & Auxiliary Water Connector Kit.

Instructions for Use

	Use-by date		Batch or lot code
	Manufacturer		Date of Manufacture
	Catalog number Product Code		Order number
	Sterilized using Ethylene Oxide		Do not re-use / Single use
	Single sterile barrier system		Product not made with natural rubber latex
	Single sterile barrier system with protective packaging inside		Consult Instructions for Use
	Do not re-sterilize		Keep dry
	Does not contain DEHP		Packing Quantity
	Do not use if package is damaged		Federal law (U.S.A.) restricts this device to sale, distribution and use by or on the order of a physician.
	Keep away from sunlight		

IMPORTANT INFORMATION

Caution: Federal law restricts this device to sale by or on the order of a physician. Read all instructions carefully before use. They contain essential information on using this device safely and effectively. Keep these instructions in a safe, accessible location, as you may need to refer to them again. If you have any questions or comments about any information in these instructions, please contact Micro-Tech USA.

INTENDED USE

The single use **Air/Water Valve (A)** is used to control the air / water function of an endoscope during GI endoscopic procedures.

The single use **Suction Valve (B)** is used to control the suction function of an endoscope during GI endoscopic procedures.

The single use **Biopsy Valve (D)** is used to cover the opening to the biopsy/suction channel of Olympus® gastrointestinal endoscopes. The Biopsy Valve provides access for endoscopic device passage and exchange, helps maintain insufflation, minimizes leakage of biomaterial from the biopsy port through the endoscopic procedure and provides access for irrigation.

The single use **Auxiliary Water Connector (C)** is used in conjunction with irrigation tubing, intended to provide irrigation via irrigation fluids such as sterile water supplied to the Olympus® GI endoscope during gastrointestinal endoscopic procedures when used in conjunction with an irrigation pump. The Auxiliary Water Connector is manufactured with a one-way valve to minimize the risk of cross-contamination of the irrigation system.

The **Air/Water Channel Cleaning Adapter** is intended as part of the precleaning process, to help clear the air/water channel of Olympus® GI endoscopes post procedures.

WARNINGS AND PRECAUTIONS

These devices are intended to be used with Olympus® GI Endoscopes **EXCLUDING Olympus® echoendoscopes that utilize a balloon channel**

- Follow the Olympus® guidelines for air/water/suction use in the Olympus® endoscopes.
- The manufacturer did not design these devices to be reprocessed or reused, and therefore cannot verify that reprocessing can clean and/or sterilize or maintain the structural integrity of these devices to ensure patient and/or user safety.
- These disposable medical devices are not intended for reuse. Any institution, practitioner, or third party, who reprocesses, refurbishes, remanufactures, reesterilizes, and or reuses these disposable medical devices must bear full responsibility for their safety and effectiveness.
- Endoscopic procedures should only be performed by persons having adequate training and familiarity with endoscopic techniques.
- Consult the medical literature relative to techniques, contraindications, complications and hazards prior to any endoscopic procedure.
- When handling these devices use appropriate aseptic technique.

Air/Water & Suction Valves:

- Make sure the holes in the air/water and suction valves are not clogged or blocked. (This could cause continuous air/CO2 flow, or inadequate suction.)
- These devices should be discarded after being used on one patient. It should not be used on multiple patients.
- The air / water valve in the endoscopic system provides backflow prevention function to the air / water channel. Not using the air / water valve can cause potential contamination to the air / water system.

Biopsy Valve:

- Do not use a sharp or pointed object to prime the valve prior to use.

- Do not leave a device hanging from the valve. Doing so can cause the creation of a larger valve slit/hole that may cause leakage.
- If the lid of the valve is opened while attached to the endoscope during a procedure, scope suction will be compromised and leakage may occur.
- If leakage occurs, a sterile gauze should be used to cover the valve.
- The device should be discarded after being used on one patient. It should not be used on multiple patients.

Auxiliary Water Connector:

- Never attach or reattach the Auxiliary Water Connector to an endoscope that is/was inside the patient and has not been reprocessed.
- Do not insert this device into the patient. The device should be attached to the auxiliary water port of a flexible gastrointestinal endoscope and an irrigation source. Use proper techniques to avoid risk of aspiration of a patient.
- Always prime the auxiliary water channel or biopsy channel prior to inserting the endoscope into the patient.
- The device is part of distal irrigation system with a backflow prevention feature and is intended for single patient use only. The device should be discarded after being used on one patient. It should not be used on multiple patients.

Air/Water Channel Cleaning Adapter:

- Do Not Use the device during patient procedures. Continuous air flow could cause patient injury. In addition, insufficient aspiration and leakage cause cross-contamination.
- Endoscopic cleaning procedures should only be performed by persons having adequate training with endoscope cleaning techniques.

Product Name:

AW01-OLY4
Disposable Air/Water, Suction, Biopsy Valve & Auxiliary Water Connector Kit.

AW01-OLY5
Disposable Air/Water Channel Cleaning Adapter and Air/Water, Suction, Biopsy Valve & Auxiliary Water Connector Kit.

Production Date: See packaging

Shelf Life: 3 years

Contents supplied sterile.

This device is not made with natural rubber latex.

AW01-OLY4

Air/Water, Suction, Biopsy Valve & Auxiliary Water Connector Kit for Olympus® 160/180/190/240/260/290 Series GI Endoscopes

AW01-OLY5

Air/Water Channel Cleaning Adapter and Air/Water, Suction, Biopsy Valve & Auxiliary Water Connector Kit for Olympus® 160/180/190/240/260/290 Series GI Endoscopes

PREPARATION

1. Prior to clinical use, familiarize yourself with the device and read all the Instructions for Use.
2. Inspect the package for shipping or handling damage. If damage is evident, **DO NOT USE THESE DEVICES. SAVE THEM AND RETURN TO MANUFACTURER AND CONTACT YOUR LOCAL PRODUCT SPECIALIST.**

INSTRUCTIONS FOR USE

Suction Valve:

3. Remove the suction valve from the tray.
4. Insert the suction valve into suction cylinder of an Olympus® endoscope (see Fig. 1).
5. Twist valve slightly to seat in proper orientation within the cylinder.
6. Confirm that the valve fits securely and that there is no gap between the valve end cap and the endoscope.
Note: The valve should not be able to be rotated once it is fully inserted into the cylinder.
7. With the thumb fully on the suction valve button, depress fully downward and release as required.

Air/Water Valve:

8. Remove the air/water valve from the tray.
9. Insert the air/water valve into the air/water cylinder of an Olympus® endoscope (see Fig. 1).
10. Confirm that the valve fits securely and that there is no gap between the valve end cap and the endoscope.
11. Depress and release the valve to check if the valve moves smoothly and freely returns to its OFF position.

Biopsy Valve:

12. Remove the biopsy valve from packaging.
13. Securely place the biopsy valve onto the biopsy/suction channel opening of an Olympus® endoscope (See Fig. 1).
Note: This biopsy valve does not accommodate a luer-lock syringe without an irrigating adaptor.
14. Close the lid when there is not a need to pass a device through the channel.
15. When the device diameter is over 3.2mm, open the lid and pass the device through the hole into the biopsy/suction channel.

Auxiliary Water Connector:

16. For use of the irrigation tubing, reference instructions for use of irrigation tubing.
17. Attach the auxiliary water connector to the auxiliary water port by screwing it onto the port (See Fig. 2). When attaching the auxiliary water connector to the port, take care to ensure that the auxiliary water connector is aligned to the port and avoid excessive tightening. Misalignment or excessive tightening of the auxiliary water connector could lead to the possibility of leakage. Once the auxiliary water connector is attached to the auxiliary water port of the endoscope, attach Irrigation Tubing to the auxiliary water connector.
18. In the case there is any leakage due to misalignment or excessive tightening, please replace with a new auxiliary water connector.
19. Discard the auxiliary water connector after each patient use.




Fig. 1

Fig. 2

Air/Water Channel Cleaning Adapter: (Applicable to AW01-OLY5)

20. Turn off the airflow regulator on the light source after initial rinsing of scope (per facility guidelines).
21. Remove the air/water valve from the endoscope.
22. Insert the air/water channel cleaning adapter into the air/water cylinder on the control body of the endoscope.
23. Place the distal end of the scope into a basin of water.
24. On the light source, turn on the airflow regulator and set it to a maximum output of "high".
25. Depress the air/water channel cleaning adapter for 30 seconds.
Note: This flushes the air/water channel with water.
26. Release the air/water channel cleaning adapter.
27. Keep the distal end of the scope in the basin of water.
28. Flush air for 10 seconds.
Note: Air bubbles will be seen in the basin of water.
29. Turn off the air/flow regulator on the light source.
30. Remove distal end of endoscope from basin of water.
31. Remove the air/water channel cleaning adapter from endoscope after precleaning step is completed.

PRODUCT DISPOSAL

After use, this product may be a potential biohazard. Handle and dispose of the product in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

WARNING

An issued or revision date for these instructions is included for the user's information. In the event that two years have elapsed between this date and product use, the user should contact Micro-Tech USA to determine if additional information is available.

Report any serious incident that has occurred in relation to the device to the distributor or manufacturer or competent authority of the member state.

Olympus® is a registered trademark of Olympus Corporation.

LIMITED WARRANTY AND DISCLAIMERS

Limited Warranty to Buyer. Micro-Tech USA warrants to Buyer that, for the earlier of one (1) year from the date of purchase, or until the product is used by Buyer. The products will be free from defects in materials and workmanship when stored and used in accordance with the instructions for storage and use provided by Micro-Tech USA and in accordance with applicable regulatory requirements.

Descriptions or specifications appearing in Micro-Tech USA's literature are meant to generally describe the products and do not constitute any express warranties.

In the event that Micro-Tech USA gives technical advice with respect to the product, it is agreed that such advice is given without any liability on Micro-Tech USA's part. Any guarantee of specific properties of or in the products shall only be effective if and to the extent specifically confirmed by Micro-Tech USA in writing. These warranties shall not apply for product failure or deficiency due to improper storage, alteration, or the consequences of uses for which the products were not designed or that adversely affect the products' integrity, reliability, or performance.

Disclaimer and Release. THE WARRANTIES, OBLIGATIONS, AND LIABILITIES OF Micro-Tech USA AS SET FORTH HEREIN ARE EXCLUSIVE. BUYER HEREBY WAIVES ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, ARISING BY LAW OR OTHERWISE, WITH RESPECT TO THE PRODUCTS AND ANY OTHER GOODS OR SERVICES DELIVERED BY BUYER, INCLUDING, BUT NOT LIMITED TO: (1) ALL OTHER EXPRESS AND IMPLIED WARRANTIES, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, AND (2)

ANY IMPLIED WARRANTY ARISING FROM COURSE OF PERFORMANCE, COURSE OF DEALING, OR USAGE OF TRADE.

Implied Warranties. The purchase of products may be subject to laws in the territories applicable to the sale of the products by Micro-Tech USA to Buyer that impose implied warranties, conditions, or obligations that cannot be excluded, restricted, or modified, or can only be excluded, restricted, or modified to a limited extent. The provisions of Paragraphs 2 and 4 shall apply to the greatest extent allowed by such laws.

Limitation of Liability. Except to the extent prohibited by applicable law, Micro-Tech USA's liability under this warranty is limited to: (a) the replacement of the products or the resupply of equivalent products; (b) the repair of the products or payment of the cost of repairing the products; or (c) payment of the cost of replacing the products or acquiring equivalent products. Micro-Tech USA shall have no obligation or liability, whether arising in contract (including warranty), tort (including active, passive, or imputed negligence, strict liability, or product liability) or otherwise, for any special, consequential, punitive, incidental, or indirect damages, or for loss of use, loss of revenue, loss of business, lost profit, or other financial loss arising out of or in connection with any product or other goods or services furnished by Micro-Tech USA, even if Micro-Tech USA was aware of the possibility of such damages or loss.

Manufactured for:

Micro-Tech USA Inc.
2855 Boardwalk Drive,
Ann Arbor, MI 48104 USA
Tel: 734-259-3768
Toll free: 877-552-4027
info@micro-tech-usa.com
Made in P.R. China