

ABC™ 3-STEP INFLATION BALLOONS

For a simple, affordable solution to balloon dilation, choose the ABC Dilation Balloon. Manufactured from a high-quality non-compliant Nylon material, the ABC Balloon offers the radial strength required to accurately and consistently dilate strictures.

The balloon can be inflated in 3 Steps, eliminating the need for multiple balloons in the same procedure. It offers an economical solution versus typical radial expansion balloons.

The balloon is available as a flexible fixed wire esophageal 8cm balloon on a 6Fr catheter, featuring rapid inflation/deflation performance. Or as a wire-guided esophageal, pyloric, colonic 5.5cm balloon, on a dual channel 7Fr catheter with a .035' flexible tip guide-wire.

NEVER
COMPROMISE
ON QUALITY

DRAMATICALLY
IMPROVE YOUR
BOTTOM LINE

THROW AWAY
CONTRACTS
FOREVER

RELIABLE
SUPPLY
PARTNERSHIP

KEY BENEFITS

CHOICE OF MATERIAL

The durable nylon material delivers strength with accurate sizing. The non-compliant material provides the radial force required for consistent dilatation. Inflated in 3 steps, the ABC balloon provides an affordable alternative to radial expanding balloons.

VISIBILITY UNDER FLUOROSCOPY

All wire guided balloons feature two radiopaque markers at the balloon shoulder, allowing for accurate placement across the stricture. The dedicated inflation port permits the introduction of contrast to further enhance visibility if desired.

TAPERED TIP

The rounded tip aids in the smooth passage of the balloon across the stricture. The wire guided balloon is provided with a pre-loaded .035 guide-wire that has a soft, flexible tip to assist in placement in tortuous or altered anatomy.

SPECIFICATIONS

WIRE GUIDED ESOPHAGEAL, PYLORIC, COLONIC BALLOONS

Order Number	Henry Schein® Item Number	Description	Balloon Diameter (mm)	Balloon Length (cm)	Working Length (cm)	Minimum Channel Size (mm)	Package Units
BD20091	125-2689	TTS - Wire Guided Dilatation Balloon	6	5.5	240	2.8	1/Box
BD20181	125-2691	TTS - Wire Guided Dilatation Balloon	8	5.5	240	2.8	1/Box
BD20271	125-2693	TTS - Wire Guided Dilatation Balloon	10	5.5	240	2.8	1/Box
BD20361	125-2696	TTS - Wire Guided Dilatation Balloon	12	5.5	240	2.8	1/Box
BD20451	125-2698	TTS - Wire Guided Dilatation Balloon	14	5.5	240	2.8	1/Box
BD20631	125-2700	TTS - Wire Guided Dilatation Balloon	16	5.5	240	2.8	1/Box
BD20721	125-2702	TTS - Wire Guided Dilatation Balloon	18	5.5	240	2.8	1/Box

FIXED WIRE ESOPHAGEAL 3-STEP INFLATION BALLOONS

Order Number	Henry Schein® Item Number	Description	Balloon Diameter (mm)	Balloon Length (cm)	Working Length (cm)	Minimum Channel Size (mm)	Package Units
FW24075	128-5660	ABC™ Esophageal Dilatation Balloon	6	8	180	2.8	2/Box
FW24085	128-5661	ABC™ Esophageal Dilatation Balloon	8	8	180	2.8	2/Box
FW24095	128-5662	ABC™ Esophageal Dilatation Balloon	10	8	180	2.8	2/Box
FW24105	128-5663	ABC™ Esophageal Dilatation Balloon	12	8	180	2.8	2/Box
FW24125	128-5664	ABC™ Esophageal Dilatation Balloon	15	8	180	2.8	2/Box
FW24145	128-5665	ABC™ Esophageal Dilatation Balloon	18	8	180	2.8	2/Box
FW24155	128-5666	ABC™ Esophageal Dilatation Balloon	20	8	180	2.8	2/Box

CAN'T FIND IT?

Additional items may be available. Contact us if you can't find what you need.



EXCEPTIONAL
ENDOSCOPIC
SOLUTIONS.

2855 Boardwalk Drive
Ann Arbor, MI 48104
877.552.4027

MTENDOSCOPY.COM

CustomerCare@MTEndoscopy.com
MPW24000 Rev. 3