

EFFECTIVE LIGATION HEMORRHOIDS MULTIPLE BAND LIGATOR SET



Endoscopic banding has been demonstrated as an effective and simple treatment for internal hemorrhoids. Banding is readily performed in an ambulatory setting and can result in less pain and have a shorter recovery period than surgical hemorrhoidectomy.

The Micro-Tech Hemorrhoids Multiple Band Ligator Set is a cost-effective solution for the outpatient treatment of hemorrhoids during a routine colonoscopy.

KEY BENEFITS

EXCELLENT VISIBILITY

A clear ligating barrel with rear mounted bands optimizes the endoscopic view while maintaining suction.

BAND DESIGN

The square band design provides consistent shape and reliable retention. The set comes with four pre-loaded latex-free bands.

USER-FRIENDLY SYSTEM

The handle unit tie hole enables quick and easy loading of the pull wire. The alert band notifies the operator when only one band remains.

SPECIFICATIONS

Order Number	Henry Schein Item Number	Description	Minimum Working Channel (mm)	Working Length (cm)	Package Units
HL32161	131-9042	Hemorrhoids Multiple Band Ligator Set, 4 Latex Free Bands, compatible with 9.4 - 13 mm scope O.D.	2.8	220	1/Box