

## PROTECT, COLLECT, CLEAN GI ACCESSORIES

Micro-Tech Bite Blocks offer comfort and protection during upper endoscopic procedures. Available in 38Fr, 60Fr and 60Plus sizes, with or without dental rim.

Polyp Retrieval Traps offer improved efficiency, reduced handling and cost savings. The unique Quad trap features removable collection trays and a secondary filter for peace of mind.

Micro-Tech's Cleaning Brushes combine performance and compliance. Available in combination kits, the brush range features 6/10mm nylon brush heads designed for effective, gentle cleaning.



NEVER  
COMPROMISE  
ON QUALITY

DRAMATICALLY  
IMPROVE YOUR  
BOTTOM LINE

THROW AWAY  
CONTRACTS  
FOREVER

RELIABLE  
SUPPLY  
PARTNERSHIP

## KEY BENEFITS

### PROTECT

The Micro-Tech Bite Blocks are designed to protect both the endoscope and the patient. The side holes allow access for scope repositioning or suction while the soft elastic strap provides a secure, comfortable fit.

### COLLECT

Reduce handling and simplify specimen collection with an in-line polyp retrieval trap. The Quad trap features removal collection trays and a secondary filter. All traps are made from clear material.

### CLEAN

Combine efficient cleaning with scope protection and convenience. Micro-Tech cleaning brushes combine soft nylon brush heads with durable catheters that improve pushability. Available in convenient combination kits.

## SPECIFICATIONS

### BITE BLOCKS

Order Number	Henry Schein Item Number	Description	Package Units
BB60011	125-2687	Bite Block, Adult, with strap	100/Box
BB60021	125-2688	Bite Block, Pediatric, with strap	100/Box
BB60041	131-2210	Bite Block, Adult w/o dental rim. 60+Fr with strap	100/Box

### CATCHALL SCOPE TRANSPORT PAD

Order Number	Henry Schein Item Number	Description	Package Units
CP62121	131-2218	CatchAll Scope Transport Pad	25/Box

### POLYP TRAPS

Order Number	Henry Schein Item Number	Description	Package Units
PC54011	129-0676	Polyp Trap with 4-chambers	25/Box
PC54021	131-2219	Polyp trap single chamber	25/Box
PC55031	131-2214	Quad Polyp Trap with removable chambers	25/Box

### CLEANING BRUSHES

Order Number	Henry Schein Item Number	Description	Package Units
QL22011	131-2220	Single ended 6mm channel brush	50/Box
QL22021	131-2222	Double ended 6mm channel brush	50/Box
QL22031	131-2221	Single ended 6mm channel and attached 10mm valve brush	50/Box
QL22041	131-2215	Double ended 5mm - 10mm valve brush	50/Box
QL22051	131-2216	Single ended channel brush + double ended valve brush	50/Box
QL22061	131-2217	Double ended channel brush + double ended valve brush	50/Box

### CAN'T FIND IT?

Additional items may be available. Contact us if you can't find what you need.



**EXCEPTIONAL**  
ENDOSCOPIC  
SOLUTIONS.

2855 Boardwalk Drive  
Ann Arbor, MI 48104  
877.552.4027

**MTENDOSCOPY.COM**

CustomerCare@MTEndoscopy.com  
MPW60000 REV.2