

PRECISE CUT HOT SNARES

Micro-Tech offers a line of polypectomy snares which combine quality and cost-effectiveness. Each snare is made with a special braided wire that provides strong radial force, stability and precise cutting properties. The power connection is mounted independently from the finger ring slide on the distal end of the hand piece, preventing the power cord from moving when the snare is opened and closed improving handling. While holding the connection, the snare loop may be rotated by turning the handle if required.

SOFT OR FIRM LOOP

0.4mm WIRE

3x3 BRAIDED SOFT WIRE

POST SWIVELS FOR EASY CONNECTION

SNARE CAN BE ROTATED

SHORTER THROW

NEVER
COMPROMISE
ON QUALITY

DRAMATICALLY
IMPROVE YOUR
BOTTOM LINE

THROW AWAY
CONTRACTS
FOREVER

RELIABLE
SUPPLY
PARTNERSHIP

KEY BENEFITS

SHORTER THROW HANDLE

The handle throw is adjusted according to the snare size. Small and Medium snares have a shorter throw than Large or Extra Large snares. A shorter throw improves comfort for users with smaller hands.

BRAIDED WIRE

The 3x3 braided wire combines the tissue acquisition features of a firm snare with the flexibility and smooth deployment of a soft snare. Suitable for daily use in facilities looking to minimize inventory.

INNOVATIVE POWER CONNECTION

The handle features a dedicated power connection, separate from the finger ring slider. This feature permits unrestricted movement and rotation of the snare handle and may reduce procedural delays.

SPECIFICATIONS

HOT SNARES 3X3 BRAIDED

Order Number	Henry Schein Item Number	Description	Loop Diameter (mm)	Sheath Diameter (Fr)	Working Length (cm)	Minimum Channel Size (mm)	Package Units
PS51011	129-0186	Mini Braided Oval	10	7	230	2.8	10/Box
PS51021	129-0673	Small Braided Oval	15	7	230	2.8	10/Box
PS51031	129-0674	Medium Braided Oval	20	7	230	2.8	10/Box
PS51041	129-0675	Large Braided Oval	24	7	230	2.8	10/Box
PS51051	129-0187	Extra Large Braided Oval	30	7	230	2.8	10/Box

HOT SNARES FIRM

Order Number	Henry Schein Item Number	Description	Loop Diameter (mm)	Sheath Diameter (Fr)	Working Length (cm)	Minimum Channel Size (mm)	Package Units
PS51071	129-0188	Mini Oval Firm	10	7	230	2.8	10/Box
PS51081	129-0189	Small Oval Firm	15	7	230	2.8	10/Box
PS51091	129-0190	Medium Oval Firm	20	7	230	2.8	10/Box
PS51101	129-0191	Large Oval Firm	24	7	230	2.8	10/Box
PS51111	129-0192	Extra Large Oval Firm	30	7	230	2.8	10/Box
PS51121	129-0193	Jumbo Oval Firm	36	7	230	2.8	10/Box

HOT SNARES HEX FIRM

Order Number	Henry Schein Item Number	Description	Loop Diameter (mm)	Sheath Diameter (Fr)	Working Length (cm)	Minimum Channel Size (mm)	Package Units
PS52031	129-0194	Small Hex Firm	15	7	230	2.8	10/Box
PS52041	129-0195	Large Hex Firm	25	7	230	2.8	10/Box

CAN'T FIND IT?

Additional items may be available. Contact us if you can't find what you need.



EXCEPTIONAL
ENDOSCOPIC
SOLUTIONS.

2855 Boardwalk Drive
Ann Arbor, MI 48104
877.552.4027

MTENDOSCOPY.COM

CustomerCare@MTEndoscopy.com
MPW51000 Rev. 3